

# SARS CoV-2: Rapid Antigen Test

# SARS CoV-2: Rapid Antigen Test

## Kit components

1. Test Card
2. Extraction buffer vial
3. Tubes for Extraction buffer
4. Nozzle for extraction tubes
5. Nasopharyngeal swabs
6. Stand for extraction buffer tubes
7. Product insert



# SARS CoV-2: Rapid Antigen Test



# SARS CoV-2: Rapid Antigen Test



**Shelf life**  
24 months from  
date of  
manufacturing



**Storage**  
At 2-30° C



# SARS CoV-2: Rapid Antigen Test



## Testing Procedure



(1)




(2)




(3)



(4)




(5)



(6)



(7)



(8)



(9)




(10)




KOSHBIO


# SARS CoV-2: Rapid Antigen Test




**Negative**



**Positive**



**Invalid**



# SARS CoV-2: Rapid Antigen Test



# SARS CoV-2: Rapid Antigen Test



Certificate of  
**FREE SALE**





# SARS CoV-2: Rapid Antigen Test



External Evaluation	
Sensitivity	Specificity
100%	100%

ICMR Evaluation	
Sensitivity	Specificity
81%	95%





## KOSHBIO PRIVATE LIMITED

Plot No. 39, Block-I, DLF Industrial Estate, Phase-1,  
Faridabad, Haryana-121003

 +91-9810163415, +91-9910897158

 [info@koshbio.com](mailto:info@koshbio.com)  [www.koshbio.com](http://www.koshbio.com)